

## dsDNA Quantification

(Quantifluor™ Dye Systems; Promega # E2670)

1. Prepare standard curve:

Standard	Volume of Lambda dsDNA Standard	Volume of 1X TE Buffer (μL)	Final Concentration of dsDNA (ng/μL)
Blank	0	1000	0.0
A	20 μL DNA Std + 980 μL 1X TE	0	2.0
B	500 μL of Std A	500	1.0
C	500 μL of Std B	500	0.5
D	400 μL of Std C	600	0.2
E	500 μL of Std D	500	0.1
F	200 μL of Std E	800	0.02

2. In a regular 96-well plate, prepare a 1:100 dilution of each unknown sample (2.5 μL sample + 250 μL AE Buffer).
3. In a black 96-well plate, transfer 100 μL of each standard and diluted sample *in duplicate* (use multi-channel pipette and Tip One filter tips to triturate 4x and transfer diluted samples to plate).
4. Prepare Quantifluor™ dsDNA Dye working solution: (10 μL 200X dye + 1,990 μL 1X TE) and add 100 μL per well (use multi-channel pipette and Tip One filter tips).
5. Incubate for approximately 5 min at RT, protected from light.
6. Measure fluorescence using the GMB Protocol on the BioTek SynergyMx 96-well plate fluorimeter. (Ex. 504; Em. 531).
  - a. Log onto computer.
  - b. Open Gen5 1.09 software.
  - c. Create New Expt
  - d. Select GMB\_Quantifluor protocol
  - e. Maximize Protocol on left side of screen and adjust plate layout if necessary.
  - f. Under Plate at the top, select “Read”.
  - g. Click the Read button.
  - h. Save file in GMB folder.
  - i. Insert plate according to instructions on screen.
  - j. Click OK.
  - k. After the machine finishes reading the plate, and data appears on the screen, click the “Statistics” Tab. For “Data”, choose Conc x Dil.
  - l. Save as an excel spreadsheet by clicking on the excel icon button.  
(R:GM\GMB\Data\QUANTIFLUOR)
  - m. Log off and close drawer (little black button).