

DNA Isolation from Human Peripheral Blood (QIAamp Blood Maxi Kit – Spin Protocol #51194)

Note: This protocol assumes the investigator is beginning this with one full Yellow-Top (type A) BD Vacutainer of human blood (5-8.5 mL). This should yield between 100 and 300 uG of HMW gDNA.

Things to do before starting:

- Equilibrate samples to room temperature (from frozen: 90 min or can place @ 4°C the night before).
- Prepare a 70°C water bath. (30 min)
- Ensure that Buffer AW1, Buffer AW2 and Qiagen Protease have been prepared according to manufacturer's instructions.
- Thoroughly mix AL Buffer before use. If a precipitate has formed, re-dissolve by incubating at 56°C

Procedure:

- 1) Label the 50 mL conical tubes that the frozen vacutainers are in, 1 through 20. Twenty samples is the maximum number of samples the tabletop centrifuge will hold.
- 2) Either note the study ID# for each sample on the log sheet or transfer the sticky label from the tube to the log sheet.
- 3) Estimate volume of blood by comparing the blood tube to a calibrated vacutainer tube; note the volume on the log sheet.
- 4) Transfer contents of tube into the same 50 mL polypropylene conical centrifuge tube that the frozen vacutainer was stored in.
- 5) Bring volume to **9.5 mL** with 1X PBS.
- 6) Add **500 uL** Qiagen Protease, and ensure proper mixing after adding the enzyme by inverting 3X.
- 7) Add **12 mL** Buffer AL, and mix thoroughly by inverting the tube 15X, followed by vigorous shaking for at least 1 min over the sink. (Invert multiple tubes simultaneously by clamping them into a rack using another empty rack, grasping both racks and inverting them together). (**Note:** Do not add Qiagen Protease directly to Buffer AL.)
- 8) Incubate at 70°C for 10 min. (**Note:** longer incubation times will not adversely affect yield).
- 9) Add **10 mL** 100% ethanol to the sample, and mix by inverting 10X, followed by additional vigorous shaking over the sink.
- 10) Carefully decant half of the solution from step 7 onto the QIAamp Maxi column placed in a 50 mL centrifuge tube (provided), taking care not to moisten the rim. (**Note:** Do not over tighten caps; Always hold the closed column in an upright position as liquid may pass through the ventilation slots on the rims of the columns even if caps are secure.)
- 11) Centrifuge at 1850xg for 3 min. (**Note:** If the solution has not completely passed through the membrane, centrifuge again at a slightly higher speed).
- 12) Remove column, discard the filtrate into properly labeled hazardous waste container, and place the column back into the same 50 mL centrifuge tube. Decant the remainder of the

GMB002A

- solution from step 7 onto the column. Close cap and centrifuge again at 1850xg for 3 min. (**Note:** Wipe off any spillage from the thread of the 50 mL tube before re-inserting the column to ensure the rim of the column does not get wet; Again, if the solution has not completely passed through the membrane, centrifuge again at a slightly higher speed).
- 13) Remove column, discard filtrate into properly labeled hazardous waste container, and place the column back into the same 50 mL centrifuge tube. (**Note:** Wipe off any spillage from the thread of the 50 mL tube before re-inserting the column to ensure the rim of the column does not get wet).
 - 14) Carefully, without moistening the rim, add **5 mL** Buffer AW1 to the column and centrifuge 4000 rpm (~3500xg) for 2 min. Do not discard flow through.
 - 15) Carefully, without moistening the rim, add **5 mL** Buffer AW2 to the column and centrifuge 4000 rpm (~3500xg) for 30 min.
 - 16) Place the column in a clean 50 mL tube (provided), and pour the filtrate into properly labeled hazardous waste container and discard tube into biohazard bag. (**Note:** Wipe any spillage off the column before inserting into tube).
 - 17) Pipet **1.2 mL** Buffer AE directly onto the membrane of the column and close cap. Incubate at room temperature for 5 min, and centrifuge at 4000 rpm (~3500xg) for 4 min.
 - 18) **Reload the eluate containing the DNA onto the membrane of the column. Close the cap and incubate at room temperature for 5 min, and centrifuge at 4000 rpm (~3500xg) for 10 min.** (**Note:** less than 1.2 mL will be eluted from the column, but this has no effect on DNA yield).
 - 19) Keep samples on ice or at 4°C for further analysis (Quantification and / or QC if necessary).

Other Reagents Needed

100% Ethanol

1X PBS