

Cryopreservation of Human Tissue in OCT (TCGA Project)

Note: This protocol is intended as a way to preserve in OCT tissue in a manner that allows the specimen to be used in multiple analyses without compromising structural integrity of the tissue, or the quality of the sample for biochemical and/or genetic investigations and the subsequent tissue sectioning on a cryostat.

Things you will need:

- OCT Embedding Media
- Embedding molds, large
- 2L Dewar $\frac{3}{4}$ full with liquid nitrogen
- Freezing basket, pre-cooled in liquid nitrogen
- Disposable forceps
- Container of dry ice
- Heavy-duty aluminum foil squares, 3" x 3"
- Sample labels

- 1) Fill the bottom of the embedding mold with just enough OCT to cover the entire area. Place at -80°C or in the vapor phase of liquid nitrogen – store at -80°C until use.
- 2) Prepare a good number of squares of heavy-duty aluminum foil roughly 3" x 3".
- 3) Have the sample labels ready and pre-made (if applicable).
- 4) Before collection, be sure to have a Dewar (2L) of liquid nitrogen ready, as well as a Styrofoam container of dry ice.
- 5) Place the appropriate number of OCT molds on the dry ice and proceed to collection.
- 6) Collect tissue/biopsy specimens per approved protocol.
- 7) Keep all specimens on ice until embedded to avoid degradation.
- 8) Quickly make note of the specimen weight – TCGA requests at least 100mg of tissue.*
- 9) Using disposable forceps, carefully remove an OCT raft from the dry ice, and quickly place a single specimen directly onto the frozen OCT. Make sure the specimen is laying flat and that the largest surface is in full contact with the OCT layer.
- 10) Quickly cover the specimen with additional OCT. Use just enough OCT to completely cover the sample.
- 11) Quickly place the OCT-embedded specimen in the basket and suspend the sample in the vapors of the liquid nitrogen until frozen.
- 12) Carefully remove the frozen specimen from the basket and wrap in one of the foil squares.
- 13) Quickly apply the completed label to the wrapped specimen and place in the liquid nitrogen, or on dry ice.
- 14) Repeat as necessary. Bank all tissue.
- 15) When all of the specimens have been frozen, they should be stored at -80°C .

* If there is plenty of tissue (more than 100mg) excise a small portion to be flash-frozen (see GMB)