

FACScan

LASER	PMT#	FILTER	EM RANGE	FLUOROCHROMES
BLUE (Ar488)	FL1	530/15	522-537nm	FITC, Calcein, CFSE, Alexa488, EGFP
	FL2	585/21	575-589nm	PE, PI
	FL3	650LP	>650nm	PerCP, 7AAD, PI, PE-Cy5

Please consult the Flow Core or <http://www.bdbiosciences.com/spectra> to determine which fluorochromes cannot be used together due to high spectral overlaps.