

Melody

| LASER | PMT# | DICHROIC FILTER | EM RANGE |
|-----------------------|------|-----------------|-----------|
| BLUE (Ar488) | A | 665LP 700/54 | 673-727nm |
| | B | 507LP 527/32 | 511-543nm |
| | C | | |
| RED (640 HeNe) | A | 752LP 783/56 | 755-811nm |
| | B | - 660/10 | 655-665nm |
| GREEN (532) | A | 752LP 783/56 | 755-811nm |
| | B | 665LP 697/58 | 668-726nm |
| | C | 605LP 613/18 | 604-622nm |
| | D | 582/15 | 575-589nm |

Please consult the Flow Core or <http://www.bdbiosciences.com/spectra> to determine which fluorochromes cannot be used together due to high spectral overlaps.